

Product : Ethyl Alcohol 200 PF ACS
Catalog# : CH5036
CAS : 64-17-5
Formula : C₂H₅OH
MW : 46.07
Grade : ACS

TEST	SPECIFICATION
Appearance	Clear colorless liquid
Assay	≥ 99.5%
Water	≤ 0.2%
Color (APHA)	≤ 10
Solubility in Water	≤ To Pass Test
Residue on Evaporation	≤ 0.001%
Acetone/Isopropyl Alcohol	≤ To Pass Tests
Titration Acid	≤ 0.0005 meq/g
Titration Base	≤ 0.0002 meq/g
Methanol	≤ 0.1%
Substances Darkened by Sulfuric Acid	To Pass Test
Substances Reducing Permanganate	To Pass Test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

