

Product : Ethylene Glycol Reagent  
Catalog# : CH5039  
CAS : 107-21-1  
Formula : C<sub>2</sub>H<sub>6</sub>O<sub>2</sub>  
MW : 62.07  
Grade : Reagent

TEST	SPECIFICATION
Appearance	Colorless to Faint Yellow liquid
Infrared Spectrum	Conforms to Structure
Purity (GC)	≥ 99.0 %
Color (APHA)	≤ 10
Acidity	≤ 0.005 %
Water (H <sub>2</sub> O)	≤ 0.3 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

