

Product : Hexane  
Catalogue : CH5046  
CAS : 100-54-3  
Formula :  $\text{CH}_3(\text{CH}_2)_4\text{CH}_3$   
MW : 86.18 g/mol  
Grade : ACS

Test	Specification
Appearance	Colorless liquid.
Density (g/ml at 25°C)	0.653 to 0.673
Color (APHA)	≤ 10 APHA
Purity (GC)	≥ 98.5%
Residue on evaporation	≤ 0.0003%
Sulfur compounds (as S)	≤ 0.005%
Thiophene	To pass the test
Water (by Karl fisher)	≤ 0.001%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

