

Product : Isopropyl Alcohol Lab  
Catalog# : CH5051  
CAS : 67-63-0  
Formula : (CH<sub>3</sub>)<sub>2</sub>CHOH  
MW : 60.10

TEST	SPECIFICATION
Appearance (Color)	Colorless liquid
Color (APHA)	≤ 10 APHA
Purity (GC)	≥ 99.0 %
Residue on Evaporation	≤ 9 ppm
Carbonyl compounds	≤ 0.002%
Acidity	To pass the test
Boiling range	82.3 – 83.3 °C
Solubility in water	To pass the test
Titration acid or base	NMT 0.001 meq/g
Water (by Karl Fischer)	≤ 0.2 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

