

Product : Isopropyl alcohol
 Catalogue : CH5053
 CAS : 67-63-0
 Formula : (CH₃)₂CHOH
 MW : 60.10
 Grade : ACS

Test	Specification
Appearance (Color)	Colorless liquid
Color (APHA)	10 APHA
Purity (GC)	≥ 99.5 %
Residue on Evaporation	≤ 10 ppm
Specific Gravity at 25°C	0.7830 to 0.7840
Carbonyl compounds	≤ 0.002%
Acidity	To pass the test
Boiling range	82.3 – 83.3 °C
Solubility in water	Pass ACS test
Titration acid or base	NMT 0.0001 meq/g
Water (by Karl Fischer)	≤ 0.2 %
Nickel (Ni)	≤ 0.1 ppm
Copper (Cu)	≤ 0.1 ppm
Iron (Fe)	≤ 0.1 ppm
Heavy metal (as Pb)	≤ 0.2 ppm

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

