

Product : Isopropyl ether ACS
Catalog# : CH5054
CAS : 108-20-3
Formula : (C₃H₇)₂O
MW : 102.18
Grade : ACS

TEST	SPECIFICATION
Appearance	Clear liquid
Assay	≥ 99.0 %
Color (APHA)	≤ 25
Peroxide (as C ₂ H ₁₄ O ₂)	≤ 0.05 %
Residue after evaporation	≤ 0.01 %
Titration acid	≤ 0.0007 meq/g

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

