

Product : Methanol
Catalogue : CH5055
CAS : 67-56-1
Formula : CH₃OH
MW : 32.04
Grade : ACS

Test	Specification
Appearance	Colorless liquid.
Color (IN APHA)	≤ 10 APHA
Purity (By GC)	≥ 99.8%
Water	≤ 0.1 %
Residue on evaporation	≤ 0.001 %
Titration acid	≤ 0.0003 meq/g
Titration base	≤ 0.0002 meq/g
Carbonyl compounds: Acetone	≤ 0.001%
Carbonyl compounds: Formaldehyde	≤ 0.001%
Carbonyl compounds: Acetaldehyde	≤ 0.001%
Substances reducing permanganate	To Pass the test
Substances darkened by sulfuric acid	To Pass the test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

