

Product : 2-Methoxyethanol ACS (Ethylene glycol monomethyl ester)
Catalog# : CH5056
CAS : 109-86-4
Formula : CH₃OCH₂CH₂OH
MW : 76.10
Grade : ACS

TEST	SPECIFICATION
Appearance	Clear liquid
Assay	≥ 99.3 %
Color (APHA)	≤ 10
Titration acid	≤ 0.002 meq/g
Water (H ₂ O)	≤ 0.1 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

