

Product : Methyl Alcohol Anhydrous ACS
Catalog# : CH5057
CAS : 67-56-1
Formula : CH₃OH
MW : 32.04
Grade : ACS

TEST	SPECIFICATION
Appearance	Clear, colorless liquid
Assay	≥ 99.8 %
Substances darkened by sulfuric acid	Passes test
Substances reducing permanganate	Passes test
Solubility in water	Passes test
Color (APHA)	≤ 10
Water (H ₂ O)	≤ 0.1 %
Residue after evaporation	≤ 0.001 %
Carbonyl	≤ 0.001 %
Titration acid	0.0003 meq/g
Titration base	00002 meq/g

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

