

Product : 1-Methyl-2-Pyrrolidinone Biotech
Catalog# : CH5062
CAS : 872-50-4
Formula : C₅H₉NO
MW : 99.13
Grade : Biotech

TEST	SPECIFICATION
Appearance	Colorless liquid
Infrared Spectrum	Conforms to Structure
UV Absorbance 400nm	≤ 0.01
UV Absorbance 350nm	≤ 0.01
UV Absorbance 325nm	≤ 0.10
UV Absorbance 300nm	≤ 0.50
UV Absorbance 285nm	≤ 1.0
Purity (GC)	≥ 99.70 %
Water (H ₂ O)	≤ 0.005 %
Color (APHA)	≤ 20
Free Amines % as Methylamine	≤ 0.01
Titrateable Base	≤ 0.001

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

