

Product : 1-Methyl-2-Pyrrolidinone ACS
Catalog# : CH5063
CAS : 872-50-4
Formula : C₅H₉NO
MW : 99.13
Grade : ACS

TEST	SPECIFICATION
Appearance	Clear Liquid
Assay	≥ 99.0%
Color (APHA)	≤ 50
Chloride content:	≤ 1 ppm
Impurity content	≤ 0.01 %
Water	≤ 0.05 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

