

Product : Methyl Tert-butyl Ether (MTBE)
Catalog# : CH5064
CAS : 1634-04-4
Formula : CH₃OC(CH₃)₃
MW : 88.15

TEST	SPECIFICATION
Appearance	Clear liquid
Assay	≥ 99.0 %
Color (APHA)	≤ 10
Peroxide (as H ₂ O ₂)	≤ 1 ppm
Residue after evaporation	≤ 0.001%
Water (H ₂ O)	≤ 0.05 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

