

Product : n-Octyl Alcohol Lab (1-Octanol)  
Catalog# : CH5067  
CAS : 111-87-5  
Formula : C<sub>8</sub>H<sub>18</sub>O  
MW : 130.23

TEST	SPECIFICATION
Appearance	Clear colorless liquid; penetrating odor
Assay	≥ 98.5 %
Color (APHA)	≤ 9
Residue after evaporation	≤ 0.04 %
Titration acid	≤ 0.002 meq/g

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

