

Product : n-Propyl Alcohol ACS  
Catalog# : CH5069  
CAS : 71-23-8  
Formula : CH<sub>3</sub>(CH<sub>2</sub>)<sub>2</sub>OH  
MW : 60.10  
Grade : ACS

TEST	SPECIFICATION
Appearance	
Assay	≥ 99.5%
Solubility in water	≤ To pass the test
Color (APHA)	≤ 10
Residue after evaporation	≤ 0.001%
Titration acid	≤ 0.0004 meq/g
Ethyl Alcohol	≤ 0.01%
Methanol	≤ 0.01%
Isopropyl Alcohol	≤ 0.05%
Carbonyl Compounds	≤ 0.03%
Water	≤ 0.2%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

