

Product : Pentane  
Catalogue : CH5070  
CAS : 109-66-0  
Formula :  $\text{CH}_3(\text{CH}_2)_3\text{CH}_3$   
MW : 72.15 g/mol  
Grade : -

Test	Specification
Appearance	Colorless Liquid
Infrared spectrum	Conforms to Structure
Refractive index at 20°C	1.356 to 1.360
Purity (GC)	≥ 99.0 %
Residue on Evaporation	≤ 0.0005 %
Water (By Karl Fischer)	≤ 0.02%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

