

Product : 1-Pentanol ACS
Catalog# : CH5071
CAS : 71-41-0
Formula : $\text{CH}_3(\text{CH}_2)_3\text{CH}_2\text{OH}$
MW : 88.15
Grade : ACS

TEST	SPECIFICATION
Appearance	liquid, with mild characteristic odor
Aasay	$\geq 98.0 \%$
Water	$\leq 0.5 \%$
Color (APHA)	≤ 30
Residue after evaporation	$\leq 0.003 \%$
Acid and esters	$\leq 0.075 \text{ meq/g}$
Carbonyl compounds (as HCHO)	≤ 0.1

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

