

Product : Petroleum ether
Catalogue : CH5072
CAS : 8032-32-4
Formula : -
MW : 82.2
Grade : ACS

Test	Specification
Appearance	Clear liquid with characteristic order.
Acidity	To Pass Test
Color (APHA)	10 max. APHA
Doctor test	To Pass Test
Final boiling temperature at 760 mm Hg	Max. 60 °C
Initial boiling temperature at 760 mm Hg	Min. 35 °C
Residue after Evaporation	Max. 0.001 %
Total Sulfur	Max. 5 ppm

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

