

Product : Propylene glycol
Catalogue : CH5075
CAS : 57-55-6
Formula : C₃H₈O₂
MW : 76.09
Grade : ACS

Test	Specification
Appearance	Colorless liquid.
Color (IN APHA)	≤ 10 APHA
Purity (By GC)	≥ 99.5%
Water	≤ 0.20 %
Residue on ignition	≤ 0.005 %
Titration acid	≤ 0.0005 meq/g
Chloride	≤ 0.0001%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

