

Product : Tetrahydrofuran lab (THF)
Catalog# : CH5078
CAS : 109-99-9
Formula : C₄H₈O
MW : 72.11

TEST	SPECIFICATION
Appearance	Clear, colorless liquid
Assay	≥ 98.5 %
Color (APHA)	≤ 19
Peroxide (as H ₂ O ₂)	≤ 0.015 %
Residue after evaporation	≤ 0.03 %
Water (H ₂ O)	≤ 0.05 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

