

Product : Toluene
Catalogue : CH5080
CAS : 108-88-3
Formula : C₆H₅CH₃
MW : 92.14
Grade : ACS

Test	Specification
Appearance	Clear, Colorless Liquid
Assay	NLT 99.5 %
Color (APHA)	10 Max APHA
Residue after Evaporation	0.001% Max
Sub. Darkened by Sulfuric Acid (H ₂ SO ₄)	To Pass Test
Sulfur compounds (as S)	0.003% max
Water (H ₂ O)	0.03% max

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

