

Product : Toluene Lab
Catalog# : CH5081
CAS : 108-88-3
Formula : C₆H₅CH₃
MW : 92.14

TEST	SPECIFICATION
Appearance	Clear, Colorless Liquid
Assay	≥ 99.0 %
Color (APHA)	≤ 10 APHA
Residue after Evaporation	≤ 0.01%
Sulfur compounds (as S)	≤ 0.003%
Water (H ₂ O)	≤ 0.03%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

