

Product : Xylenes
Catalogue : CH5083
CAS : 1330-20-7
Formula : C₆H₄(CH₃)₂
MW : 106.17 g/mol
Grade : -

Test	Specification
Appearance	Clear colorless liquid.
Assay (Xylene Isomers & Ethyl benzene) corrected for water	NLT 98.5%
Color (APHA)	10 max
Substances darkened by sulfuric acid	To pass the test
Residue after evaporation	≤ 0.002%
Sulfur compounds (as S)	0.003 % max
Water (By KF)	≤ 0.05%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

