

Product : Acetic Acid, Glacial, Lab
Catalog# : CH8002
CAS : 64-19-7
Formula : CH₃COOH
MW : 60.05

TEST	SPECIFICATION
Appearance	Liquid with a sharp odor
Assay	≥ 99.0 %
Color (APHA)	≤ 10
Dilution test	Passes test
Residue after evaporation	≤ 0.01%
Acetic anhydride [(CH ₃ CO) ₂ O]	≤ 0.001%
Chloride (Cl)	≤ 2ppm
Sulfate (SO ₄)	≤ 1ppm
Substance reducing dichromate	≤ passes test
Substance reducing permanganate	≤ passes test
Titration base	≤ 0.004 meq/g

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

