

Product : Acetic Acid. Glacial, FCC
Catalog# : CH8003
CAS : 64-19-7
Formula : CH₃COOH
MW : 60.05
Grade : FCC

TEST	SPECIFICATION
Assay (C ₂ H ₄ O ₂)	99.5-100.5 wt.%
Lead	≤ 0.5 mg/kg
Nonvolatile Residue	≤ 0.005%
Readily Oxidizable Substances	To Pass Test
Solidification Point	≥ 15.6°C

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

