

Product : Ammonium Hydroxide Lab  
Catalog# : CH8007  
CAS : 1336-21-6  
Formula : NH<sub>4</sub>OH  
MW : 35.05  
Grade : Lab

TEST	SPECIFICATION
Appearance	Clear liquid with a suffocating odor
Assay (as NH <sub>3</sub> )	27.0 – 30.0 %
Residue after ignition	≤ 0.02%
Carbon Dioxide (CO <sub>2</sub> )	≤ 0.02%
Chloride (Cl)	0.4 ppm
Nitrate (NO <sub>3</sub> )	1 ppm
Phosphate (PO <sub>4</sub> )	1 ppm
Sulphate (SO <sub>4</sub> )	1 ppm
Substances reducing permanganate	Passes test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

