

Product : Ammonium Hydroxide, ACS, Trace Grade
 Catalog# : CH8009
 CAS : 1336-21-6
 Formula : NH₄OH
 MW : 35.05
 Grade : ACS

TEST	SPECIFICATION
Appearance	Clear, colorless liquid with a suffocating odor
Assay (as NH ₃)	20-30 %
Chloride (Cl)	0.5 ppm
Phosphate (PO ₄)	2 ppm
Sulfate (SO ₄)	2 ppm
Mercury (Hg)	1 ppb
Aluminum (Al)	1 ppb
Barium (Ba)	1 ppb
Boron (B)	5 ppb
Cadmium (Cd)	1 ppb
Calcium (Ca)	1 ppb
Chromium (Cr)	1 ppb
Cobalt (Co)	1 ppb
Copper (Cu)	1 ppb
Iron (Fe)	5 ppb
Lead (Pb)	1 ppb
Lithium (Li)	1 ppb
Magnesium (Mg)	1 ppb
Manganese (Mn)	1 ppb
Molybdenum (Mo)	1 ppb
Potassium (K)	1 ppb
Silicon (Si)	5 ppb
Sodium (Na)	5 ppb
Strontium (Sr)	1 ppb
Tin (Sn)	1 ppb
Titanium (Ti)	1 ppb



Product : Ammonium Hydroxide, ACS, Trace Grade
Catalog# : CH8009
CAS : 1336-21-6
Formula : NH₄OH
MW : 35.05
Grade : ACS

Test	Specification
Vanadium (V)	1 ppb
Zinc (Zn)	1 ppb
Zirconium (Zr)	1 ppb

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

