

Product : Formaldehyde Solution ACS
Catalog# : CH8015
CAS : 50-00-0
Formula : HCHO
MW : 30.03
Grade : ACS

TEST	SPECIFICATION
Appearance	Colorless solution
Assay	36.5-38.0%
Color (APHA)	≤ 10
Residue after ignition	≤ 0.005%
Titration acid	≤ 0.006meq/g
Chloride (Cl)	≤ 5ppm
Sulfate (SO ₄)	≤ 0.002%
Heavy metal (by ICP-OES)	≤ 5ppm
Iron (Fe)	≤ 5ppm

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

