

Product : Formic acid 88% Lab  
Catalog# : CH8018  
CAS : 64-18-6  
Formula : HCOOH  
MW : 46.03  
Grade : Lab

TEST	SPECIFICATION
Appearance	Clear, colorless liquid
Assay	≥ 88.0%
Color (APHA)	≤ 15
Dilution test	Passes test
Residue after evaporation	≤ 0.02%
Acetic acid (CH <sub>3</sub> COOH)	≤ 0.4%
Ammonium (NH <sub>4</sub> )	≤ 0.05%
Chloride (Cl)	≤ 0.01%
Sulfate (SO <sub>4</sub> )	≤ 0.02%
Sulfite (SO <sub>3</sub> )	Passes test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

