

Product : Hydriodic acid 55% ACS
Catalog# : CH8021
CAS : 10034-85-2
Formula : HI
MW : 127.91
Grade : ACS

TEST	SPECIFICATION
Appearance	Clear liquid; may darken with time
Assay	≥ 55.0- 58.0%
Free iodine (I ₂)	≤ 0.75%
Residue after ignition	≤ 0.01%
Chloride and bromide (as Cl)	≤ 0.05%
Sulfate (SO ₄)	≤ 0.005%
Phosphate (PO ₄)	≤ 0.001%
Heavy metals (as Pb)	≤ 0.001%
Iron (Fe)	≤ 0.001%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

