

Product : Hydrobromic Acid Reagent
Catalog# : CH8023
CAS : 10035-10-6
Formula : HBr
MW : 80.91

TEST	SPECIFICATION
Appearance	Clear liquid, colorless to pale yellow when fresh; may become yellow to brown on storage
Assay	46.5-48.0 %
Residue after ignition	≤ 0.02%
Chloride (Cl)	≤ 0.5%
Iodide (I)	≤ 0.03%
Phosphate (PO ₄)	≤ 0.01%
Sulfate and sulfite (as SO ₄)	≤ 0.03%
Selenium (Se)	≤ 0.1 ppm

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

