

Product : Hydrochloric Acid 30% Lab
Catalog# : CH8024
CAS : 7647-01-0
Formula : HCl
MW : 36.46
Grade : Lab

TEST	SPECIFICATION
Appearance	Clear fuming liquid
Assay	29.5% - 31.0%
Color (APHA)	≤ 10
Residue after ignition	≤ 5 ppm
Sulfate (SO ₄)	≤ 1 ppm
Free Chloride (Cl)	≤ 1 ppm
Sulfite (SO ₃)	≤ 1 ppm

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

