

Product : Hydrochloric Acid 37% ACS
Catalog# : CH8025
CAS : 7647-01-0
Formula : HCL
MW : 36.46
Grade : ACS

TEST	SPECIFICATION
Appearance	Clear, fuming liquid
Assay	36.5-38.0 %
Color (APHA)	≤ 10
Residue after ignition	≤ 5ppm
Bromide (Br)	≤ 0.005 %
Sulfate (SO ₄)	≤ 1ppm
Sulfite (SO ₃)	≤ 1ppm
Free chlorine (Cl)	≤ 1ppm
Ammonium (NH ₄)	≤ 3ppm
Arsenic (As)	≤ 0.01 ppm
Heavy metals	≤ 1ppm
Iron (Fe)	≤ 0.2 ppm

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

