

Product : Hydrofluoric Acid 48% ACS
Catalog# : CH8026
CAS : 7664-39-3
Formula : HF
MW : 20.01
Grade : ACS

TEST	SPECIFICATION
Appearance	Clear liquid
Assay	48.0-51.0 %
Fluosilicic acid (H ₂ SiF ₆)	≤ 0.01%
Residue after ignition	≤ 5 ppm
Chloride (Cl)	≤ 5 ppm
Phosphate (PO ₄)	≤ 1 ppm
Sulfate and sulfite (as SO ₄)	≤ 5 ppm
Arsenic (As)	≤ 0.05 ppm
Copper (Cu)	≤ 0.1 ppm
Iron (Fe)	≤ 1 ppm
Heavy metals (as Pb)	≤ 0.5 ppm

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

