

Product : Hydrofluoric Acid 70% Lab
Catalog# : CH8029
CAS : 7664-39-3
Formula : HF
MW : 20.01

TEST	SPECIFICATION
Appearance	Clear liquid
Assay	68.0% - 71.0 %
Fluosilicic acid (H ₂ SiF ₆)	≤ 0.01%
Residue after ignition	≤ 4 ppm
Chloride (Cl)	≤ 4 ppm
Phosphate (PO ₄)	≤ 1 ppm
Sulfate and sulfite (as SO ₄)	≤ 4 ppm

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

