

Product : Hydrogen peroxide 30% ACS
Catalog# : CH8030
CAS : 7722-84-1
Formula : H₂O₂
MW : 34.01
Grade : ACS

NO₃ %

TEST	SPECIFICATION
Appearance	Clear liquid
Assay	29.0-32.0 %
Color (APHA)	≤ 10
Residue after evaporation	≤ 0.002 %
Titration acid	≤ 0.0006 meq/g
Chloride (Cl)	≤ 2ppm
Nitrate (NO ₃)	≤ 2ppm
Phosphate (PO ₄)	≤ 2ppm
Sulfate (SO ₄)	≤ 5ppm
Ammonium (NH ₄)	≤ 5ppm
Heavy metals (as Pb)	≤ 1ppm
Iron (Fe)	≤ 0.5ppm

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

