

Product : Hydrogen peroxide 6%
Catalog# : CH8036
CAS : 7722-84-1
Formula : H₂O₂
MW : 34.01

TEST	SPECIFICATION
Appearance	Clear liquid
Assay	5.5% - 7.0%
Color (APHA)	≤ 10
Residue after evaporation	≤ 0.02%
Chloride (Cl)	≤ 2 ppm
Phosphate (PO ₄)	≤ 2 ppm
Sulfate (SO ₄)	≤ 4 ppm

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

