

Product : Molybdic Acid ACS
Catalog# : CH8038
CAS : 7782-91-4
Formula : -
MW : -
Grade : ACS

TEST	SPECIFICATION
Appearance	White or slightly yellowish solid
Assay	≥ 85.0%
Insoluble in dilute ammonium hydroxide	≤ 0.01%
Chloride (Cl)	≤ 0.002%
Arsenate, phosphate and silicate (as SiO ₂)	≤ 0.001%
Phosphate (PO ₄)	5 ppm
Sulfate (SO ₄)	≤ 0.2%
Heavy metals (as Pb)	≤ 0.003%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

