

Product : Nitric acid 70% ACS  
Catalog# : CH8040  
CAS : 7697-37-2  
Formula : HNO<sub>3</sub>  
MW : 63.01  
Grade : ACS

TEST	SPECIFICATION
Appearance	Clear colorless liquid, fuming in moist air; with a characteristic odor
Assay	68.0-70.0 %
Color (APHA)	≤ 10
Residue after ignition	≤ 5ppm
Chloride (Cl)	≤ 0.5ppm
Sulfate (SO <sub>4</sub> )	≤ 1ppm
Arsenic (As)	≤ 0.01ppm
Heavy metals (as Pb)	≤ 0.2ppm
Iron (Fe)	≤ 0.2ppm

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

