

Product : Nitric Acid 90% Fuming, ACS
Catalog# : CH8041
CAS : 7697-37-2
Formula : HNO₃
MW : 63.01
Grade : ACS

TEST	SPECIFICATION
Appearance	Clear, yellow fuming liquid
Assay	≥ 90.0%
Dilution test	Passes test
Residue after ignition	≤ 0.002%
Dissolved oxide	≤ 0.1 % asN ₂ O ₃
Chloride (Cl)	0.7 ppm
Sulfate (SO ₄)	5 ppm
Arsenic (As)	0.3 ppm
Heavy metals (as Pb)	5 ppm
Iron (Fe)	2 ppm

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

