

Product : Perchloric Acid 60% ACS
Catalog# : CH8043
CAS : 7601-90-3
Formula : HClO₄
MW : 100.46
Grade : ACS

TEST	SPECIFICATION
Appearance	Clear, colorless liquid
Assay	60.0-62.0%
Color (APHA)	≤ 10
Residue after ignition	≤ 0.003%
Silicate and phosphate (as SiO ₂)	5 ppm
Chloride (Cl)	≤ 0.001%
Nitrogen compound (as N)	≤ 0.001%
Sulfate (SO ₄)	≤ 0.001%
Heavy metal (as Pb)	1 ppm
Iron (Fe)	1 ppm

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

