

Product : meta-Phosphoric acid ACS
Catalog# : CH8045
CAS : 37267-86-0
Formula : -
MW : -
Grade : ACS

TEST	SPECIFICATION
Appearance	Transparent, glass- like solid or very silky masses; hygroscopic
Assay	33.5-36.5 %
Stabilizer (NaPO ₃)	57.0-63.0 %
Nitrate (NO ₃)	≤ 0.001 %
Sulfate (SO ₄)	≤ 0.005 %
Chloride (Cl)	≤ 0.001 %
Arsenic (As)	≤ 1ppm
Heavy metals (as Pb)	≤ 0.005 %
Iron (Fe)	≤ 0.005 %
Substance residue permanganate (as H ₃ PO ₃)	≤ 0.02 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

