

Product : Phosphoric Acid 70% Lab
Catalog# : CH8046
CAS : 7664-38-2
Formula : H_3PO_4
MW : 98.0
Grade :

TEST	SPECIFICATION
Appearance	Clear, viscous liquid
Assay	≥ 70.0 %
Color (APHA)	≤ 10
Insoluble matter	≤ 0.01 %
Chloride (Cl)	≤ 3ppm
Nitrate (NO ₃)	≤ 5ppm
Sulfate (SO ₄)	≤ 0.003 %
Potassium (K)	≤ 0.05 %
Sodium (Na)	≤ 0.025 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

