

Product : Phosphoric Acid 85% Food Grade
Catalog# : CH8049
CAS : 7664-38-2
Formula : H₃PO₄
MW : 98.00

TEST	SPECIFICATION
Assay (H ₃ PO ₄) wt.%	85.0% - 85.5%
Specific gravity @ 25/15.5°C	1.686 - 1.692
Arsenic, as As, ppm	≤ 3 max
Cadmium, as Cd, ppm	≤ 3 max
Chloride, as Cl, ppm	≤ 10 max
Color, APHA	≤ 10 max
Fluoride, as F, ppm	≤ 10 max
Iron as Fe, ppm	≤ 10 max
Lead, as Pb, ppm	≤ 3 max
Sulfate, as SO ₄ , ppm	≤ 250 max

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

