

Product : Phosphoric Acid 85% Lab
Catalog# : CH8050
CAS : 7664-38-2
Formula : H₃PO₄
MW : 98.00

TEST	SPECIFICATION
Appearance	Clear, viscous liquid
Assay	≥ 85.0%
Color (APHA)	≤ 10
Insoluble matter	≤ 0.01%
Chloride (Cl)	≤ 2 ppm
Sulfate (SO ₄)	≤ 0.03%
Reducing substance	Passes test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

