

Product : Phosphotungstic Acid Crystal Reagent
Catalog# : CH8051
CAS : 12501-23-4
Formula : $H_5O_{41}PW_{12}$
MW : 2880.05
Grade : Reagent

TEST	SPECIFICATION
Appearance	White to Off-White Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity) 1 g + 10 mL, H ₂ O	Clear to Slightly Hazy
Loss on Drying	≤ 15 %
Sulfate (SO ₄)	≤ 0.02%
Chloride (Cl)	≤ 0.03%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

