

Product : Sulfuric Acid ACS
Catalog# : CH8053
CAS : 7664-93-9
Formula : H₂SO₄
MW : 98.07
Grade : ACS

| TEST | SPECIFICATION |
|--------------------------------------------------------|-----------------------------------------|
| Appearance | Free from suspended or insoluble matter |
| Assay | 95.0-98.0 % |
| Color (APHA) | ≥ 10 |
| Residue after ignition | ≤ 5ppm |
| Chloride (Cl) | ≤ 0.2ppm |
| Nitrate (NO ₃) | ≤ 0.5ppm |
| Ammonium (NH ₄) | ≤ 2ppm |
| Substances reducing permanganate (as SO ₂) | ≤ 2ppm |
| Arsenic (As) | ≤ 0.01ppm |
| Heavy metals (as Pb) | ≤ 1ppm |
| Iron (Fe) | ≤ 0.2ppm |
| Mercury (Hg) | ≤ 5ppb |

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

