

Product : Sulfuric Acid Lab
Catalog# : CH8054
CAS : 7664-93-9
Formula : H₂SO₄
MW : 98.07

Sulfuric Acid Lab

TEST	SPECIFICATION
Appearance	Clear, colorless liquid
Assay	95.0-98.0%
Color (APHA)	≤ 10
Residue after ignition	4 ppm
Chloride (Cl)	0.2 ppm
Substance reducing permanganate (as SO ₂)	2 ppm

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

