

Product : Trichloroacetic Acid 85%
Catalog# : CH8058
CAS : 76-03-9
Formula : CCl_3COOH
MW : 163.38

| TEST | SPECIFICATION |
|-------------------------------------|--------------------------|
| Appearance | Colorless to white solid |
| Assay | $\geq 85.0\%$ |
| Clarity of solution | Passes test |
| Insoluble matter | $\leq 0.1\%$ |
| Residue after ignition | $\leq 0.3\%$ |
| Chloride (Cl) | $\leq 0.02\%$ |
| Phosphate (PO_4) | 4 ppm |
| Sulfate (SO_4) | $\leq 0.02\%$ |
| Substance darkened by sulfuric acid | Passes test |

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

