

Product : Trichloroacetic Acid ACS  
Catalog# : CH8060  
CAS : 76-03-9  
Formula : CCl<sub>3</sub>COOH  
MW : 163.38  
Grade : ACS

TEST	SPECIFICATION
Appearance	Colorless to white solid
Assay	≥ 99 .0%
Clarity of solution	Passes test
Insoluble matter	≤ 0.01%
Residue after ignition	≤ 0.03%
Chloride (Cl)	≤ 0.002%
Nitrate (NO <sub>3</sub> )	≤ 0.002%
Phosphate (PO <sub>4</sub> )	≤ 5 ppm
Sulfate (SO <sub>4</sub> )	≤ 0.02%
Heavy metals (as Pb)	≤ 0.002%
Iron (Fe)	≤ 0.001%
Substance darkened by sulfuric acid	Passes test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

