

Product : Hydrogen peroxide 30% Lab
Catalog# : CH8064
CAS : 772-84-1
Formula : 34.01
MW : H₂O₂

TEST	SPECIFICATION
Appearance	Clear liquid
Assay	25.0-35.0 %
Chloride (Cl)	≤ 3ppm
Nitrate (NO ₃)	≤ 2ppm
Phosphate (PO ₄)	≤ 2ppm
Sulfate (SO ₄)	≤ 5ppm
Calcium (Ca)	≤ 1ppb

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

